

SRD

SEDNA Dispersing Impeller



• Features



고속의 RPM에서 사용되며,
파우더를 분산, 용해하는 것에 적합합니다.

Running at high RPM,
it is suitable for dispersing and dissolving powder.

적용분야 시럽, 점안제, 현탁액, 파우더 등

Applications Syrup, Eye drop, Suspension, Powder etc.

• Product Information

	Model	Mixing volume Design RPM	Specification
	SRD-60	50L (AC 1.5kw) 1500 RPM	<ul style="list-style-type: none"> ▶ Impeller <ul style="list-style-type: none"> • Material : Stainless Steel, Hastelloy , Titanium • Surface Finish : Ra <0.4μm
	SRD-100	150L (AC 2.2kw) 1500 RPM	<ul style="list-style-type: none"> ▶ Welding Plate <ul style="list-style-type: none"> • Material : Stainless Steel, Hastelloy, Titanium • Surface Finish : Ra <0.4μm ▶ Bearing <ul style="list-style-type: none"> • Bearing Material : SiC / TC
	SRD-150	700L (AC 4.0kw) 1500 RPM	<ul style="list-style-type: none"> • Elastomer : EPDM, Viton, Silicon PTFE Capsule etc. ▶ Drive Unit <ul style="list-style-type: none"> • 220/380~440 VAC • Protection : IP54/55 • 50/60Hz
	SRD-200	1500L (AC 7.5kw) 1500 RPM	
	SRD-220	3000L (AC 11kw) 1500 RPM	